Area of application for the ZESTRON® Bath Analyzer 30:



In order to achieve consistent cleaning results, the user has to ensure that the concentration of the cleaning agent stays within the recommended application range. Various factors such as "drag-out," dilution and evaporation can potentially influence the actual concentration. Regular measurements of all process related parameters are therefore highly recommended. The ZESTRON® Bath Analyzer 30 is a simple, easy to use method that provides reliable and accurate results for fresh as well as contaminated cleaning baths.

Simple measurement in only a few steps:

① Sampling / Decanting



2 Adding test solution



3 Shaking / Evaluating



Advantages of the ZESTRON® Bath Analyzer 30:

- Very precise concentration measurement for fresh as well as contaminated cleaning baths
- Simple and fast procedure
- Requires little space with minimal investment

Specifically developed for the following cleaning agents:

HYDRON® WS 400

The ZESTRON® Bath Analyzer 30 Kit includes:



- 1) PE bottle with injection nozzle
- 2) 2 pairs of vinyl gloves
- 3) Digital thermometer (°C/°F)
- 4) Sampling beaker
- 5) Glass cylinder with marking
- 6) 6 x 100 mL ZESTRON® Bath Analyzer 30 test solution
- 7) User manual, including correlation charts (not illustrated)

To place an order:

| Part No. | Product Description | |
|-----------|-------------------------------------|--|
| 0425-0000 | ZESTRON® Bath Analyzer 30 | |
| | Portable case | |
| 0426-0100 | ZESTRON® Bath Analyzer 30 | |
| | 6 x 100 mL refill bottles in carton | |
| 0426-0005 | ZESTRON® Bath Analyzer 30 | |
| | 5 L refill bottle | |

Contact: ZESTRON America 11285 Assett Loop Manassas, Virginia 20109 Phone: (703) 393-9880

Fax: (703) 393-8618 Email: infousa@zestron.com